

OXYGEN CONCENTRATOR



OLV-B1 portable oxygen concentrator with *Battery*
is the most lightest, smallest and quietest oxygen
concentrator for travel and ambulatory oxygen needs.

Olive CARE FOR YOU!

Zhengzhou Olive Electronic Technology Co.Ltd
Tel: 0086 371 8608 1306 Fax: 0086 371 6906 7368

www.zzolive.com
E-mail: info@zzolive.com

Features

1. Continuous flow, 0-5L/min adjustable.
2. Light weight: 5kg
3. Compact size: 9.6" * 7.1" * 14.3"
4. Three power suppliers: 220V/110V, 12V (Car use), Rechargeable Lithium-ion Battery.
5. Anion function
6. Fashionable carrying case with adjustable shoulder strap.
7. Low noise, less than 40db(A)



Products Specifications

Size: 9.6 x 7.1 x 14.3(Inches) or 24.5 x 18.1 x 36.5(cm)

Weight: machine 5kg, battery 1kg

Continuous Flow Settings: 0-5LPM

Purity Oxygen: $\geq 90\% \pm 3\%$ (1LPM)

Output Pressure: 40Kpa-60Kpa

Sound Level: less than 40DB(A)

Power: AC Power 220V \pm 22V, 50Hz \pm 1Hz/110V \pm 15V, 60Hz \pm 1Hz

DC Power 12V

Rechargeable Battery 7800mAh

Power Consumption: 90W

Battery Duration: 60mins

Battery recharge time: 4 hours recharge time to achieve 90% capacity

